

01

Syringes

Advantages

- > Certified to ISO 9001:2000, ISO 13485:2003 and CE mark standards.
- > Sterile single use disposable.
- > Latex free.
- > Individually packaged, simple to use.
- > We offer non-sterile large quantity in bulk case for Re-serialise Via EtO or Autoclave.
- > We also offer a different type syringe with attached needle.

3 Part Luer Lock Syringes - Without Needles

Wego Code	SPE.	PCS/CTN	Package
0110102	2ml	1600	Blister
0110103	3ml	1600	Blister
0110105	5ml	1200	Blister
0110110	10ml	800	Blister
0110120	20ml	600	Blister
0110130	30ml	400	Blister
0110150	50ml	300	Blister
0110160	60ml	300	Blister
0110100	100ml	160	Blister



3 Part Luer Slip Lock Syringes - Without Needles

Wego Code	SPE.	PCS/CTN	Package
0111101	1ml	2000	Blister
0111102	2ml	1600	Blister
0111103	3ml	1600	Blister
0111105	5ml	1200	Blister
0111110	10ml	800	Blister
0111120	20ml	600	Blister
0111130	30ml	400	Blister
0111150	50ml	300	Blister
0111160	60ml	300	Blister
0111100	100ml	160	Blister



3 Part Eccentric Tip Syringes

Wego Code	SPE.	PCS/CTN	Package
0112105	5ml	1200	Blister
0112110	10ml	800	Blister
0112120	20ml	600	Blister
0112130	30ml	400	Blister
0112150	50ml	300	Blister
0112160	60ml	300	Blister



Luer Lock Safety Syringes - With Needles

Wego Code	SPE.	PCS/CTN	Package
0113103	3ml	1000	Blister
0113105	5ml	1000	Blister
0113110	10ml	1000	Blister



Insulin/Fixed Needle Syringes

Wego Code	SPE.	PCS/CTN	Package
0114101	1ml	3000	Blister

