

02

Needles

Advantages

- > Certified to ISO 9001:2000, ISO 13485:2003 and CE mark standards.
- > Sterile single use disposable.
- > Latex free.

- > Individually packaged, simple to use.
- > We offer non-sterile large quantity in bulk case for Re-sterilise Via EtO or Autoclave.
- > We also offer a different type syringe with attached needle.

Hypodermic Needles

Wego Code	Colour	SPE.	PCS/CTN	Package
0210116	White	16G	5000	Blister
0210118	Pink	18G	5000	Blister
0210119	Yellow	19G	5000	Blister
0210120	Yellow	20G	5000	Blister
0210121	Green	21G	5000	Blister
0210122	Black	22G	5000	Blister
0210123	Blue	23G	5000	Blister
0210124	Purple	24G	5000	Blister
0210125	Orange	25G	5000	Blister
0210126	Light Blue	26G	5000	Blister
0210127	Grey	27G	5000	Blister



Blood Collection Needle

Wego Code	Colour	SPE.	PCS/CTN	Package
0211120	Yellow	20G	5000	Blister
0211121	Green	21G	5000	Blister
0211122	Black	22G	5000	Blister

Features

- > Plastic protector for intravenous needle.
- > Colour-coded caps and hubs which allow for rapid identification of needle gauge size.
- > Used together with the needle holder.
- > Scope for use; for collecting intravenous blood sample in clinic. Connect with Vacuum Blood Tube.

