

# 03

## Blood Bags

### Blood Bags - SINGLE

Wego Code	Anti-Coagulant	Capacity	Sets/Carton
0310100	CPDA-1	100ml	100
0310250	CPDA-1	250ml	100
0310350	CPDA-1	350ml	100
0310450	CPDA-1	450ml	100
0310500	CPDA-1	500ml	100



### Blood Bags - DOUBLE

Wego Code	Anti-Coagulant	Capacity	Sets/Carton
0320250	CPDA-1	250ml	80
0320350	CPDA-1	350ml	80
0320450	CPDA-1	450ml	80
0320500	CPDA-1	500ml	80



### Blood Bags - TRIPLE

Wego Code	Anti-Coagulant	Capacity	Sets/Carton
0330250	CPDA-1	250ml	50
0330350	CPDA-1	350ml	50
0330450	CPDA-1	450ml	50
0330500	CPDA-1	500ml	50
0331250	CPDA-SAGM	250ml	45
0331350	CPDA-SAGM	350ml	45
0331450	CPDA-SAGM	450ml	45



### Blood Bags - QUADRUPLE

Wego Code	Anti-Coagulant	Capacity	Sets/Carton
0340250	CPDA-1	250ml	40
0340350	CPDA-1	350ml	40
0340450	CPDA-1	450ml	40
0340500	CPDA-1	500ml	40

# 04

## Infusion & Transfusion Set

### Description

- > Used for the transfusion of blood or blood products.
- > Injection moulded and blown chamber.
- > Single packed with poly-bag or blister package.
- > Pyrogen Free.
- > Kink Free PVC Tube.
- > EO Sterilised Single Use.

### Infusion Set

Wego Code	Description
0410101	With Luer Lock
0411101	With Luer Slip
0412101	With Air-vent



### Transfusion Set

Wego Code	Description
0413101	Without Air-vent
0412102	With Air-vent & Y injection port

