

06

Scalp Vein Set

Description

- > Used for blood collection, with a sterile closed system for drawing blood from difficult veins, such as the back of the hand.
- > Butterfly shaped wings for better handling and attachment to skin.
- > Silicone cannula for smooth penetration.

- > Slip or Luer set end matches well with different types of infusion sets.
- > Kink Free with Ultra-sharp needles.
- > Tube length 30cm or other if required.
- > EO sterilised Single Use.

Scalp Vein Set			
Wego Code	Colour	SPE.	PCS/CTN
0610119	Yellow	19G	2000
0610121	Teal	21G	2000
0610122	Brown	22G	2000
0610123	Blue	23G	2000
0610125	Orange	25G	2000



07

I.V. Cannula

Features

- > Certified to ISO 9001:2000, ISO 13485:2003 and CE mark standards.
- > Advanced tip designed protection cap with a recessed plug to cover the injection port.
- > Stainless steel silicone tipped needle for smooth penetration.
- > Thin wall construction permits higher flow rates.
- > EO Sterilised Single Use.
- > Latex Free, Non-Pyrogenic.
- > Individual blister packaging, simple to use.

With Wings			
Wego Code	Colour	SPE.	PCS/CTN
0711114	Yellow	14G	200
0711116	Grey	16G	200
0711118	Teal	18G	200
0711120	Purple	20G	200
0711122	Dark Blue	22G	200
0711124	Yellow	24G	200



Y or I Type			
Wego Code	Colour	SPE.	PCS/CTN
0710118	Blue	18G	200
0710120	Pink	20G	200
0710122	Dark Blue	22G	200
0710124	Yellow	24G	200

